




SVED™ Wound Treatment System



Benefits of the SVED™ Wound Treatment System:

- Negative Pressure Wound Therapy to promote wound healing
 - Provides NPWT for removal of fluids and exudate
 - Helps maintain a moist wound environment
 - Helps keep wounds clean and removes infectious materials
- Continuous or intermittent pressure
- Open cell hydrophilic foam dressing
- Wound irrigation feature provides therapy options and patient comfort
- SpeedConnect™ for ease of application
- One - touch pressure adjustment with lock outs
- Visual and audible dressing leak alarms
- Automatic shutoff prevents canister overflow

The SVED™ System can be used with various aqueous topical solutions*

* When applying topical solutions, be sure to apply the product consistent with the solution manufacturer's instructions for use.

Ordering Information

Part #	Description
6701132	SVED™ Therapy Unit (weight: 2 lbs.)
47-4000	SVED™ 300cc Canister Set with Gel (10 per case)
47-1702	Svamp™ Small Dressing Set (10 per case)
47-1701	Svamp™ Medium Dressing Set (10 per case)
47-1700	Svamp™ Large Dressing Set (10 per case)
47-1703	Svamp™ XL Dressing Set (10 per case)
47-1751	Vit™ Foam Dressing (10 per case)
47-2000	SpeedConnect™ Tubing Set (10 per case)
47-1760	ACS™ - Abdominal Dressing Set (5 per case)
47-2500	Y-Connector (10 per bag)
46-307B	Drape 10" x 13" (10 per case)
Also Available:	
6700032	SVEDMAN™ Therapy Unit (weight: 6 lbs.)
71-3ITI-4	SVEDMAN™ 1200cc Canister Set with Gel (10 per case)

Dressing set complete with Svamp™ foam dressing, tubing, connectors, drape and irrigation set.

DISTRIBUTED BY:



691 Green Crest Drive
 Columbus, OH 43081
 phone: 866-MED-BEDS (633-2337)
 fax: 614-237-1177
www.amemedbeds.com